



PL-6018II

Trolley Color Doppler

- Display: 15" LCD monitor.
- Gray scale: 256.
- Probe frequency: 2.0-12Mhz.
- Mode: B, BB, 4B, left & right B/M, M, PW, B/D, CFM, B/C, BC/B, B/BC.
- Duplex: real-time simultaneous 2D,
- Doppler, frequency simultaneously in three modes cape real-time status.
- Cine-loop: 1000 frame cine loop memory.
- Storage media: 64G massive image-storage capacity.

Focusing:

- Electronic focusing + acoustic lens focusing.
- Transmission focus: single focus, double focus, 3-focus and 4- focus.
- Continuously dynamic focusing mode is used for front end reception.



Software Package:

Basic measurement package, obstetrics measurement package, andrology measurement package, gynecology measurement package, urology measurement package, vascular measurement package, small parts package, orthopedic measurement package, cardiac measurement package.

Measurement and Calculations:

- B-mode basic measurements: distance, angle, perimeter and area (ellipse method, track method), volume histograms, cross section.
- M-mode basic measurements: heart rate, time, distance, speed.
- Doppler measurement: time, heart rate, speed, acceleration.

Memory Function :

- Probe parameters memory
- Image memory
- Cine memory
- Measurement result memory
- Report memory

Others:

- USB ports: video, S-video, VGA, 4 USB ports.
- Probe connector: 2 active transducer connector ports.
- Power supply: (AC)220V±22V 50Hz±1Hz.

Standard Configuration:

- Scanner Body: 1pc
- 3.5Mhz Convex Probe: 1pc
- Charger with a Set of Wires: 1set
- Bottle of USG gel: 1pc
- User Manual: 1pc

Optional:

- 5.0Mhz Micro-convex Probe
- 6.5MHz Trans-vaginal Probe
- 7.5Mhz Linear Probe



Xuzhou Palmary Electronics Co., Ltd



ADD.: No.4 Building, Zhongkai City Light, West Huaihai Road, Xuzhou City, Jiangsu Province, China (221008)
 TEL.: 86-516-80221766 / 83892858 FAX: 86-516-83892808 / 83892898
 BUSINESS: ultrasoundscanner@palmary.cn SERVICE: paierdianzi@palmary.cn
 www.palmaryultrasound.com